Lobe Rotor Pump BCL



The BCL series present lobe rotor pumps designed according to the EHEDG specifications that comply with the strictest hygienic requirements of the biotechnology and pharmaceutical industry.

The applications include pumping of viscous as well as low-viscous media that require a gentle transfer. These pumps are also used in the applications requiring perfect cleanness and sterilisation by CIP/SIP processes.

The aseptic version is especially indispensable for the applications where a sterile pumped product is to be safe from any contamination.

## Operating principle

The SLR pumpss basically consist of two lobe rotors which rotate synchronously inside a casing without touching each other.

As the rotors rotate, the spaces between the lobes and the casing are successively filled with fluid, which is transported to the discharge nozzle with a fixed amount of displacement. The pumped fluid forms a continuous stream thanks to the tolerances between the lobes and the pump casing, thus ensuring an efficient pumping.

The operating principle is exactly the same as the one of the SLR lobe rotor pump.

Moreover, the pump is completely drainable as it is mounted vertically.

Apart from a double mechanical seal aseptic design pumps have steam barriers at the front cover as well as at the clamp connections of the inlet and outlet of the pump.

Thus all places where the product may have contact with the atmosphere have barriers.

### Design and features

Vertical support. Completely drainable pump.

Clamp connections.

Tri-lobe rotors.

Interior surface roughness Ra<0.5  $\mu m$ .

EPDM gaskets according to FDA 177.2600 normative.

Hygienic single mechanical seal.

Parts in contact with the media: AISI 316L.

The assembly of the lobes is effected wihitout screws in contact with the media.



#### Materials

Parts in contact with the media:

Support:

Gaskets (standard)

Mechanical seal (standard)

AISI 316L GG-15

EPDM in compliance with FDA

TYPE EL

SiC/SiC/EPDM







## **Options**

Single mechanical seal SiC/C/EPDM Single cooled mechanical seal (quench) Back to back mechanical seal FPM and PTFE gaskets

Aseptic version:

steam barrier at the pump cover steam barriers at the pump inlet and outlet connections double mechanical seal SiC/SiC/EPDM-SiC/C/EPDM

Traceability and Certification of the materials EN 10204 3.1


Surface finish certificate


ATEX design available


Stainless steel support (CF-8)

## Technical specifications

Max. flow40 m3/hMax. head10 barMax. working temperature120 °CMax. speed950 rpm







BCL	DN	da	db	dc	de	ea	eb	ec	ed	ee	sp	tb	vb	VC	vd	ve	vf	vh	za	zb	zc
1-25	1"	165	112	25	160	40	6	21,6	19	6	71	287	115	135	9	85	10	105	31	162	225
1-40	1 1/2"										77	299									231
2-40	1 1/2"	200	140	31	190	50	8	27	24	8	76	339	125	150	11	105	12	130	51,5	177	263
2-50	2"										82	351									269
3-50	2"	280	190	46,5	250	80	10	41,4	38	12	91	430	170	210	13	130	14	170	74	232	348
3-80	3"	200									101	452									360

## Насосы — Роторный насос BCL

# http://k-tep.com.ua/

# ⊠ k-tep@ukr.net

■ Office +38 044 2091823 ■ MTC +38 066 9076563 ■ Киевстар +38 098 3676414

Skype: k-teppumps

### Замечания

Условия поставки: DDP склад г. Киев

### Заметки

Время поставки рассчитано согласно дате предложения и изменяется в зависимости от даты подтверждения заказа.

Помещая заказ, Покупатель принимает предложение и все спецификации, характеристики и условия, указанные в данном документе.