Tank bottom aseptic valve

Sanitary
flow
equipment

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

NET VOLUME OF VALVE BODY CAVITY WITH PTFE DIAPHRAGM INSTALLED

Tank bottom valve body available in the following type:


VFA10°-S/TC

10° Angled butt weld


VFA10°-TK/TC

Angled TK removable body


VFA- S/TC

Flanged flush weld

VFA- FL/TC

Flanged removable body


Net volume in ml

| Valve Code/Size | INCH | DN | VFA10°-S/TC | VFA10°-TK/TC | VFA- S/TC | VFA- FL/TC |
|------------------|--------|-----|--|--------------|-----------|---------------|
| VFA10°-S/TC1" | 1" | 25 | 14 | 50 | 25 | 30 |
| VFA10°-S/TC1"1/2 | 1" 1/2 | 40 | 27 | 72 | 75 | 105 |
| VFA10°-S/TC2" | 2" | 50 | 92 | 235 | 85 | 150 |
| VFA10°-S/TC3" | 3" | 80 | 170 | 310 | 320 | Not available |
| VFA10°-S/TC4" | 4" | 100 | Not standard product, available on specific request only | | | |

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
Special executions available on request:

- welding plate with radius or bigger different thickness, special adaptors plate
- 45°Outlet port butt weld
- Stainless steel handle/pneumatic actuator
- Adjustable flow regulator
- Mechanical/inductive control box

AVAILABLE RADIAL DIAPHRAGMS

STANDARD

Material: SILICONE
Code: MVA-X
Available size: 3/4"-2"

Material: TFM 1705 PTFE
Code: MVA-G
Available size: 3/4"-4"

BELLOW ON REQUEST

Material: TFM 1705 PTFE
Code: MSVA-G

Double stroke for high viscosity media

Material: TFM+INOX
Code: MSVA-G INOX

Aisi 316L safety coating cap with unique O-Ring in FEP for save the TFM from abrasive crystal

Regulatory compliance:

- FDA 21CFR177.1550
- USP Class VI<87> and <88> (70°C and 121° C)
- ADIF animal derived ingredient free

TFM is a registered trademark of Dyneon

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The valve body are designed according to PED Directive 97/23/EC for Europe, ASME VIII Div.2 for US and F.E.M. (Finite Element Method) calculated, approved and certified by notified body.
 Warning : The applied diaphragm and actuator may have a different design temperature and/or pressure.
 The weakest part in the assembled product set the final, permitted design temperature and pressure limits.

Flow rate

In order to design valves for a process system correctly, the valve size is determined by the required flow rate. The Kv-value serves as a calculation basis for the different process conditions. This value is stated in the following table with regard to nominal diameter and standards.

Kv-value (m3/h)

The Kv-value is a parameter defining the flow rate of valves. It describes the amount of water from 5° to 25°C which flows through the valve at a pressure loss of 1 bar when the valve is 100% open


Conversion

For the correct Kv to Cv conversion calculation, use only the stated units formulas above.

The Kv-value must be converted from (cubic meter/hour) by utilizing the following conversion factor.

In the US the flow rate of water is measured with the Cv-value in US-gallons per minute (gpm) with a pressure drop of 1 PSI.

conversion of Kv to Cv $Cv = 1,17 \times Kv$
 conversion of Cv to Kv $Kv = 0,86 \times Cv$

Flow characteristics with TFM 1600 PTFE diaphragm installed (Flow In to Out)

| Valve Code/Size | INCH | DN | Kv-value (m3/h) | Cv-value (gpm) |
|-----------------|--------|-----|--|----------------|
| VFA10°-S/TC1" | 1" | 25 | 10,5 | 12,2 |
| VFA10°-S/TC1" | 1" 1/2 | 40 | 15,0 | 17,5 |
| VFA10°-S/TC2" | 2" | 50 | 32,0 | 37,4 |
| VFA10°-S/TC3" | 3" | 80 | 78,0 | 91,2 |
| VFA10°-S/TC4" | 4" | 100 | Not standard product, available on specific request only | |

Available configurations

STANDARD DESIGN
 The valves is composed of three following main parts:
 • **Actuator** - either manual or pneumatic
 • **Diaphragms** - available in Silicone and TFM
 • **Valve body** - in line, zero static, tank bottom
 The modular design ensures a wide flexibility and you can combine the three parts in almost any combination.
 The manual / pneumatic actuators and body are assembled by a locking ring

Minox
STAINLESS
STEEL
HANDLE

M
STANDARD
PLASTIC
HANDLE

RP
ADJUSTABLE
TRAVEL
SCREW

PS
MECHANICAL /
INDUCTIVE
POSITION
INDICATOR
OPEN-CLOSED

LOCKING RING

P
STAINLESS
STEEL
PNEUMATIC
ACTUATOR

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